



# HME™

## Introducing Rice Lake's WHEELCHAIR SCALE

The Rice Lake single ramp wheelchair scale allows high capacity accurate weighing for a wide range of patients. The convenient keypad tare reduces errors and increases efficiency. Rice Lake wheelchair scales come in a variety of configurations to ensure safe accurate weighing for patients needing the assistance.



large LCD display



foldable ramp

Foldable and lightweight, this sturdy scale transports easily on heavy duty rear wheels to any location.



- Two wheels for portability
- Single ramp
- Non-slip platform
- Large 1-inch LCD display
- Motiontrap™ movement technology
- Unit of measure toggle or locking
- EMR capable
- Serial and USB data ports
- Bluetooth® BLE 4.0
- Tare function
- Hold function
- BMI (body mass index) function
- Factory calibrated
- Easy to store, folded depth is 11.75 in

### Specifications

CAPACITY	DIMENSIONS	POWER	BATT TYPE	BATT USAGE	WARRANTY
1000 x 0.2 lb (450 x 0.1 kg)	29.5 x 32.0 x 2.4in	9 VDC	6 AA alkaline	25 hours	2 year

### Wheelchair Scale

MODEL#	CAPACITY	COMMUNICATION	PART#
350-10-2	1000 x 0.2 lb (450 x 0.1 kg)	USB	141445
350-10-2BLE	1000 x 0.2 lb (450 x 0.1 kg)	Bluetooth BLE 40	194739

### Dimensions

PART	INCHES	CENTIMETERS	PART	INCHES	CENTIMETERS
Scale base L	29.5	74.9	Handle base H	9.3	23.5
Scale base W	32	81.3	Scale/Handle W	37.5	95.3
Scale base H	2.4	6.1	Handle Height	45	114.3

### Options

DESCRIPTION	PART#
9 VDC adapter, 120 V/230 V auto switching c-MET* Medical grade power supply	115655
9 VDC adapter, 120 V/230 V auto switching*	133077
9 VDC adapter, 230 V**	102508
Epson TM-U220 tape printer	81076
Cable, RS232 scale to printer, 6 ft, 25-pin D-Sub (male) to RJ45	186830
Cable, USB 6 ft USB-A (male) to USB-B (male)	79074
Cable, Communication DB9 (female) to RJ45 (male), 4 conductor	217091

### Shipping Weight

MODEL#	PART#	L (in)	W (in)	H (in)	PACKAGE WGT (lbs)	DIM WGT (lbs)	SHIPPING WGT (lbs)
350-10-2	141445	42	41	14	90	174	174
350-10-2BLE	194739	42	41	14	90	174	174